


SEMPERSHIELD®

VS

PURPLE

sempermed®  
A MEMBER OF THE SEMPERIT-GROUP



Kimberly-Clark



100 GLOVES PER BOX (90 XXL)

100 GLOVES PER BOX

BOXES COLOR CODED BY SIZE

BOXES ARE NOT COLOR CODED BY SIZE

COLOR: AZURE BLUE

COLOR: PURPLE

FINGERTIP MIL THICKNESS 6.0 MILS

FINGERTIP MIL THICKNESS 5.9 MILS

PALM MIL THICKNESS 5.2 MILS

PALM MIL THICKNESS 4.7 MILS

ELONGATION: 726%

ELONGATION: 550%

TENSILE STRENGTH: 20.8 MPa

TENSILE STRENGTH: 21.0 MPa

SIZES: S - XXL

SIZES: XS - XL

NITRILE / POWDER FREE / TEXTURED

NITRILE / POWDER FREE / TEXTURED



sempermed®

A MEMBER OF THE SEMPERIT-GROUP