

Our Gloves are Disposable. The People they Protect are Not.™

SemperSure®

Accelerator Free

NITRILE EXAM

Fingertip Textured



FEATURES

- Accelerator Free
- Tested for use with chemotherapy drugs
- Chlorinated

BENEFITS

- Low dermatitis potential
- Appropriate for handling cytotoxic chemicals
- Tape won't stick to glove

REORDER INFORMATION

Glove Color: Cobalt Blue

SUNF201	X-Small	200/BX	10 BX/CS
SUNF202	Small	200/BX	10 BX/CS
SUNF203	Medium	200/BX	10 BX/CS
SUNF204	Large	200/BX	10 BX/CS
SUNF205	X-Large	180/BX	10 BX/CS
SPK-SUNF	1 Pair Each Size: Small, Medium, Large		



RECOMMENDED USES



MEDICAL



DENTAL



VET



LAB



FOOD

PRODUCT SPECIFICATIONS



LATEX
FREE



CHEMO
TESTED



TEXTURED



BEADED
CUFF



SINGLE
USE ONLY



AMBIDEXTROUS



CHLORINATED

