

Ambu[®] Man Basic



Key Features

- Robust and reliable
- Accurate training of chest compressions
- Visible compression depth feedback
- Ambu hygienic system for mouth-to-mouth and mouth-to-nose ventilation
- Used for bag-mask ventilation
- Placement of guedel airway (when using the guedel face piece)
- Realistic ventilation resistance
- Visible chest rise during ventilation
- Anatomically correct closed torso
- Easy AED training

AmbuMan Basic

The AmbuMan Basic is a cost efficient training manikin that is designed for BLS algorithm training.

The closed shaped torso (anatomically correct) allows trainees to learn and train on the use of Automatic External Defibrillators (AEDs) and the placement of pads.

The AmbuMan Basic allows for an optimal training of chest compression and ventilation. The unique hygienic system prevents cross-contamination.

The correct compression depth and required ventilation volume can be monitored by the feedback indicator on the manikin.

The AmbuMan Basic is also available as a four pack for group training.

Realistic ventilation

Trainees can feel ventilation resistance and expiration through mouth and nose. The AmbuMan Basic features look, listen, and feel feedback.

Cost efficient and reliable

The AmbuMan Basic is made to be cost efficient and reliable for BLS training.

Class training

For easy transportation and set up the AmbuMan Basic comes in a soft bag either packed as a single manikin or in a set of four.

Specifications

Weight Torso	6 kg
Length	70 cm
Width	37 cm
Weight 4 x torso and bag	30 kg
Max. ventilation volume	1200 ml
Max. compression depth	7.5 cm



AmbuMan Basic

Materials

Basic unit	ABS (Acrylnitril-butadien-styrol)
Head	PVC (Polyvinylchlorid) hard
Chest skin	PVC (Polyvinylchlorid) soft
Face piece	PVC (Polyvinylchlorid) soft
Headbag, disposable	PE (Polyethylene)
Training suit	Polyester
Carrying bag	PVC-coated nylon



Placement of AED pads

Order information

Item no	Description
A268 005 000	AmbuMan Basic incl. training jacket, carrying case, 1 face piece and 25 head bags
A268 006 000	AmbuMan Basic 4 pack incl. 4 AmbuMan Basic Torso with training jacket, carrying case, 4 face piece and 100 head bags
Accessories	
234 000 702	Ambu Head bags, disposable (set of 100)
234 000 703	Ambu face piece (Set of 5)
234 000 725	Ambu gudel face piece



Compression indicator