



Apex Pro™ XP100 offers excellent contrast to synthetic opiates. Law Enforcement, Corrections and EMS will benefit from this added visual protection. ApexPro™ XP100 features our 2-Ply Technology with Black Exterior and Gray Interior. Apex Pro™ XP100 features our proprietary pH balanced interior for skin health.

ApexPro™ XP100 Powder-Free Nitrile Exam Gloves were developed to offer the following:

- 2-Ply Technology with Black Exterior and Gray interior
- Black Exterior ideal for contrasting hazardous residues
- Optimized for Law Enforcement, Corrections and EMS
- 12” length, 8.25mil fingertip thickness
- Proprietary triple-wash & pH application for maintaining skin health
- Fully-textured for proven grip enhancement in dry and fluid environments
- UL certified NFPA 1999-2018 edition
- **Fentanyl Permeation Testing:** APEXPro™ XP100 has zero breakthrough with Fentanyl Citrate 100mcg/2ml for 240 minutes (test limit) per ASTM D6978. This assures Law Enforcement, Corrections and EMS maximum protection in the event of inadvertent contact with synthetic opiates. Test results on file.

Ordering information

Size	Catalog no.	Packaging
XSmall	APXP12-1X	100 gloves / box
Small	APXP12-1	100 gloves / box
Medium	APXP12-2	100 gloves / box
Large	APXP12-3	100 gloves / box
XLarge	APXP12-4	100 gloves / box
XXLarge	APXP12-5	80 gloves / box
XXXLarge	APXP12-6	80 gloves / box

All packaging – 10 boxes / case

Product specification

Specifications	ASTM D6319	NFPA 2018 Ed	ApexPro™ XP100
Length (in)	9.05”	n/a	12”
Fingertip (mils)	1.9 mil	n/a	8.25 mil
Palm (mils)	1.9 mil	n/a	6.0 mil
Cuff (mils)	n/a	n/a	4.5 mil

Tensile strength

Before aging (MPa)	14	14	22
After aging (MPa)	14	14	20

Elasticity / elongation

Before aging (MPa)	500%	500%	560%
After aging (MPa)	400%	500%	530%
Puncture resistant	n/a	1.01 lbs	1.9 lbs

