

Ambu® Baby Manikin

The hygienic lifelike baby manikin.



Key Benefits

- Patented Hygienic system
- Airway open only in sniff position
- Airway obstruction possibility
- The Ambu® Baby is a lifelike manikin that stimulates babies up to the age of one year.

Ambu® Baby Manikin

A Realistic training tool

The baby manikin corresponds to a real baby up to one year old. It is the youngest member of the Ambu family of traditional manikins. Teaching BLS with the baby manikin gives a very lifelike simulation.

Emergency training made easy

Like all Ambu BLS manikins, the baby manikin has the patented hygienic system that eliminates the risk of cross contamination. All trainees receive their own face piece and head bag for the training. At the end of the training the face piece can be cleaned and the head bag is discarded. A truly 100% hygienic way of training.

Enabling you to provide the best care

Just as in real life the airway opens only when the head is properly moved into the sniff position. The instructor has the ability to simulate an obstruction of the airway. The baby manikin has additional features such as a realistic brachial pulse, to provide correct and realistic training. Stomach ventilation can be checked visually.

Specifications

DIMENSIONS	
Length	16" (40 cm.)
Weight	5.5 lbs. (2.5 kg.)

Ordering Specifications

ITEM NO.	DESCRIPTION
256 001 000	Ambu® Baby Manikin includes 5 face pieces and 100 head bags
256 000 702	Baby head bags, 100/set
256 000 703	Baby face pieces, 5/set



Ambu head bags



Ambu Baby face pieces